




BIORAD[®]

NCM QUAT

Endoscope And Medical
Instrument Disinfection



FORMS OF PACKAGE

Description	Package
1 litre PE bottle	20 x 1 litre
5 litre PE drum	4 x 5 litre
10 litre PE drum	1 piece

CHEMICAL COMPOSITION:

Didecyldimethylammoniumchloride, Benzalkoniumchloride, corrosion and foam inhibitors and auxiliary agents.

PRODUCT SPECIFICATIONS:

NCM QUAT is an aldehyde-free and phenol-free, combined and concentrate solution activated with quaternary ammonium compounds to disinfect medical equipment. It is compliant with all materials and effective against all species of bacteria, fungi, and viruses.

INTENDED USE:

NCM Quat cleans and disinfects all types of interventional surgical equipment, dental surgical equipment, surgical laparoscopy equipment, flexible and rigid endoscopes. It is compliant with glass, ceramic, silicone, plastic, acrylic, polymer, polycarbonate, polystyrene, HDPE, porcelain and stainless steel. It is suitable for use with endoscope washing machines, ultrasonic washing machines and for manual use.

USAGE:

It is diluted by 1% for disinfection of medical instruments (e.g. in order to prepare a solution of 10L, add 100ml solution to water of 9.9L.) Action time is 5 minutes. It has maximum action time of 15 days. It is recommended to replace it if solution forms organic waste and gets dirty noticeably.

ANTI-MICROBIAL CHARACTERISTICS:

Bactericidal, Fungicidal, Virucidal, Tuberculocidal, Micobactericidal