




BIORAD®

POWDER PA

Endoscope And Medical
Instrument Disinfection



FORMS OF PACKAGE

Description	Package
HDPE bucket	1 x 10 kg
HDPE bucket	20 x 1 kg
HDPE bucket	4 x 5 kg

CHEMICAL COMPOSITION:

Tetraacetythylenediamine, Sodium Perborate, Enzymes (Amylase, Lipase, Protease) Corrosion inhibitors and auxiliary agents.

INTENDED USE:

It provides high level disinfection of all types of interventional surgical equipment, dental surgical equipment, surgical laparoscopy equipment, flexible and rigid endoscopes. It does not cause any harm on the endoscope covers such as rust, corrosion, softening and colour change. It is suitable for disinfection of medical equipment made of plastic, PVC, glass, silicone, ceramic, acrylic, polycarbonate, polymer, polystyrene, HDPE, porcelain and stainless steel.

USAGE:

It is diluted by 1% for use. Stir solution of 1% and BIORAD POWDER PA 10 g (1 level dosing spoon come together) in lukewarm of 900 ml (20-25°C) until it becomes homogenous. After activation for 10 minutes, the prepared solution is now ready for use. Prepared solution can be used for 24 hours; however, you should replace solution when it gets dirty excessively.

ANTI-MICROBIAL CHARACTERISTICS:

Bactericidal, Fungicidal, Sporocidal, Virucidal, Tuberculocidal, Micobactericidal