



# BIORAD®

## POWERFUL PA

Endoscope And Medical Instrument Disinfection



### FORMS OF PACKAGE

Description	Package
1 litre PE bottle	20 x 1 litre
5 litre PE drum	4 x 5 litre
10 litre PE drum	1 piece

### CHEMICAL COMPOSITION:

Peracetic acid, Hydrogen peroxide, acetic Acid, anticorrosive agents and auxiliary materials.

### SCOPE AND INTENDED USE:

Ensures high-level disinfection of all kinds of equipment for interventional surgery, dental surgical equipment, surgical laparoscopy equipment, flexible and rigid endoscopes. It does not cause any harm on the endoscope covers such as rust, corrosion, softening, colour change, etc. It is suitable for disinfection of medical equipment made of plastic, PVC, glass, silicone, ceramic, acrylic, polycarbonate, polymer polystyrene, HDPE, porcelain and stainless steel. It can be used with any trademarks of endoscope. And it can be washed with endoscope washing machines and it is suitable for manual use.

### USAGE:

It is diluted by 3% for use. Use dosing pump to prepare the solution. Fill water into the disinfection tank. Add Biorad Powerful Pa referring to the dilution table. Disassemble the equipment precleared with enzyme and place in the container in such a way that all surfaces and cavities will contact with solution. Seal the tank and wait for as much as the action time and then take the equipment from the solution and rinse with ample water. Efficiency of the solution should be checked with test strips. It is recommended to replace the solution when it gets dirty noticeably. Please comply with the user instructions given by the manufacturer when you use it with automatic endoscope washing and disinfection machines.

### ANTI-MICROBIAL CHARACTERISTICS:

Bactericidal, Fungicidal, Sporocidal, Virucidal, Tuberculocidal, Micobactericidal.