



FAMED SAĞLIK TIBBİ CİHAZLAR
İTH. İHRC. SAN. ve TİC. LTD. ŞTİ.



Go and Diagnose

Portable Gastroscope System



Portable Gastroscope System with Ambulatory Flight- Case

www.famed.com.tr

info@famed.com.tr

Ambulatory Flight- Case

Features

- ⦿ Robust flight case for use with TK-VETKIT gastroscopes GU Series
- ⦿ Available to hold tablet or laptop within 15 inch



LED Light Source L-10





Features

- ⦿ High brightness LED, good color restoration and long operation life.
- ⦿ Touch-screen control, friendly interactive interface.
- ⦿ Brightness with 9 levels adjustment, suitable for a variety of scenes.
- ⦿ All in one light source + water/air pump.

Specification

Model	Power Voltage	Light	Air Pump Pressure	Diameter	Brightness	Air Pump
FMD-10	100-240V AC	LED	45Kpa-65Kpa	φ8.5mm	9 levels	3 levels


Application



- ▶ Upper and lower GI tract work in companion animals (gastroscope and colonoscope) and airways investigations in equine(upper airways scope)

- ▶ Ideal for numerous large animal applications, including gastroscopy, colonoscopy, esophagoscopy, tracheoscopy, hysteroscopy, rhinoscopy, pharyngoscopy and laryngoscopy.

- ▶ Ideal for numerous large animal applications, including gastroscopy, colonoscopy, esophagoscopy, tracheoscopy, hysteroscopy, rhinoscopy, pharyngoscopy and laryngoscopy.




Specification

Model	Outer Diameter	Working Channel	Working Length	Up/Down	Left/Right	View Angle
FMD-160H	6.8mm	2.2mm	1500mm	U210°, D130°	L100°, R100°	120°
FMD-180H	8.0mm	2.2mm	1500mm	U210°, D130°	L100°, R100°	120°
FMD-190H	9.2mm	2.8mm	1500mm	U210°, D130°	L100°, R100°	120°
FMD-190HL	9.2mm	2.8mm	3000mm	U210°, D130°	L100°, R100°	120°
FMD-1B0H	11.5mm	3.2mm	3500mm	U180°, D180°	L160°, R160°	120°



Flexible & Compatible



*USB Solution, compatible with your own
PC/Laptop/ Window System Tablet*

Convenient & Reliable

- ▶ A reliable partner for veterinarians to take the gastroscop system everywhere: veterinary hospitals, farms, horse ranches, zoo, Aquarium etc.
- ▶ Make outcalls and emergency treatments possible and convenient.



FAMED SAĞLIK TIBBİ CİHAZLAR
İTH. İHRC. SAN. ve TİC. LTD. ŞTİ.



Address:: Cevizli Mh. Hacılar Cd. No:44/A Kartal - İstanbul / TURKEY
Phone : +90 216 597 0 567 / +90 216 441 29 99 **Fax:** +90 216 397 15 15
W: www.famed.com.tr **E:** info@famed.com.tr

