

FMD6D



Veterinary Anesthesia System



Friendly Powerful Reliable

Application

The anesthesia system makes a good performance in operation room with acuity from high to low, applicable for cases from simple to complex. Our expertise in anesthesia delivery is based on rich 14-year history of developing anesthesia solutions that meet your needs. The good performance is based on high security, high accuracy, high stability and accurate monitoring.



Flowmeter



Pressure gauge



APL valve



Breathing bag

Friendly Powerful Reliable

FMD6D

Veterinary System

Trust point

- The equipment could adapt to different animals and procedures freely.
- Designed and manufactured with over 14 years experience in this area.
- Flexible configurations to satisfy your needs.
- International standard and advanced technology suitable for a wide range of use.
- Over 2,000 units installed in the world.



Feature

- The high precision of flowmeter instantly lets you know how much fresh gas flow to animals.
- Integrated breathing circuit design ensures easy operating and keeps tidy.
- Vaporizer with temperature & flow compensation and self-lock function keeps the safety at anytime.


Specification

Flow	O ₂ (0 ~ 5 L/min)
Rapid oxygen supply	25 L/min ~ 75 L/min

Monitoring parameter	
Airway pressure	0 cmH ₂ O ~ 100 cmH ₂ O

Working condition	
Gas source	O ₂
Pressure	280 kPa ~ 600 kPa

Vaporizer



FAMED SAĞLIK TIBBİ CİHAZLAR İTH.İHRC. ve SAN. TİC. LTD. ŞTİ.

Address: Cevizli Mh. Hacılar Cd. No:44/A Kartal / İstanbul / TURKEY

Phone : +90 216 597 0 567 +90 216 441 29 99 Fax: +90 216 397 15 15

Web : www.famed.com.tr / Email: info@famed.com.tr