



Veterinary Disposable Sclerotherapy Needle

Precise puncture in all areas.

These disposable injection needles make sure that you get the highest reliability at injections such as in endoscopic mucosal resection. The stainless steel needle with a XX precision sharpening provides an optimal and easy puncture. The transparent PTFE—tube also enhances the rigidity, which optimizes the control of the needle once more. To prevent damage to the working channel, the needle is secured in the tube.

Offers the injection needles in different lengths and with different diameters.

Sharp needle edge
Stretched out quickly

Infuse smoothly


Glass III registration certificate Exclusive patented design

Special features

- Needle with extremely sharp precision cut
- Handle with ergonomic design
- Safety lock for controlled needle handling
- Different needle lengths and diameters
- Slif, transparent PTFE tube to prevent kinking

User-friendly control handle with safety lock.

The slim and thus particularly ergonomic handle offers maximum control of the needle. The safety—lock prevents inadvertent motion of the needle while the insertion of the instrument. Damage on the working channel will be safely avoided.



Model	Diameter of the outer tube (mm)	Length (mm)	Diameter (mm)	Length (mm)
Injection Needle				
FMD-Z5Z-18x1200x25x4	1.8	1200	25G	4
FMD-Z5Z-18x1800x25x6	1.8	1800	25G	6
FMD-Z5Z-23x1600x23x5	2.3	1600	23G	5
FMD-Z5Z-23x2300x23x5	2.3	2300	23G	5
FMD-Z5Z-23x1800x21x4	2.3	1800	21G	4
FMD-Z5Z-23x1800x21x5	2.3	1800	21G	5
FMD-Z5Z-23x1800x21x6	2.3	1800	21G	6
FMD-Z5Z-23x1800x23x4	2.3	1800	23G	4
FMD-Z5Z-23x1800x23x5	2.3	1800	23G	5
FMD-Z5Z-23x1800x23x6	2.3	1800	23G	6
FMD-Z5Z-23x1800x25x4	2.3	1800	25G	4
FMD-Z5Z-23x1800x25x5	2.3	1800	25G	5
FMD-Z5Z-23x1800x25x6	2.3	1800	25G	6
FMD-Z5Z-23x2300x21x4	2.3	2300	21G	4
FMD-Z5Z-23x2300x21x5	2.3	2300	21G	5
FMD-Z5Z-23x2300x21x6	2.3	2300	21G	6
FMD-Z5Z-23x2300x23x4	2.3	2300	23G	4
FMD-Z5Z-23x2300x23x5	2.3	2300	23G	5
FMD-Z5Z-23x2300x23x6	2.3	2300	23G	6
FMD-Z5Z-23x2300x25x4	2.3	2300	25G	4
FMD-Z5Z-23x2300x25x5	2.3	2300	25G	5
FMD-Z5Z-23x2300x25x6	2.3	2300	25G	6



Cevizli Mah. Hacilar Cad.
No:44/A Kartal / İstanbul



+90 216 597 0 567
+90 216 441 29 99



+90 216 397 15 15

www.famed.com.tr