



FAMED SAĞLIK TIBBİ CİHAZLAR
İTH. İHRC. SAN. ve TİC. LTD. ŞTİ.

FMD-2023

Veterinary Endoscopy Comprehensive Solution



**PRODUCT
CATALOG**




www.famed.com.tr



info@famed.com.tr

Veterinary Gastroscope System GU Series



21.5 inches
image processor


Gastroscope

LED light source

Trolley



Features



Work Station

- Image Processor: 21.5 inches touch screen medical grade monitor
- Work Station: Restore medical data truly and accurately
- High brightness, high contrast ratio and wide viewing angle
- Professional medical report management system

Light Source



- High brightness LED, good color restoration and long operation life
- Touch screen control, friendly interactive interface
- Brightness with 9 level adjustment, suitable for a variety of scenes
- All in one light source + water/air pump

Gastroscope



- High definition system, CMOS technology with resolution of 720P
- Ergonomic design of endoscope handle
- More efficient cooperation with 4-way angulation
- Wide angle view 120°

Advantage




Medical Record Management System

- Restore medical report truly and accurately, support report printing



Excellent Manipulation

- Ergonomics design
- Touch screen intelligent functions
- Multi-function buttons on endoscope handle



High Definition 2.0Mega pixel

- CMOS technology with 1080P resolution for better viewing experience

Technical Parameters

GU Series

Model	Outer Diameter	Working Channel	Working Length	Up/Down	Left/Right	View Angle	Function Keys	Water Jet	Resolution
FMD-160H	6.8mm	2.2mm	1500mm	U210°/D130°	L100°/R100°	120°	1	/	HD
FMD-170Q	7.8mm	2.8mm	1500mm	U180°/D130°	L100°/R100°	120°	3	/	HD
FMD-170FH	7.8mm	2.2mm	1500mm	U180°/D130°	L100°/R100°	140°	3	/	FHD
FMD-190FHP	9.2mm	2.8mm	1500mm	U180°/D130°	L100°/R100°	140°	3	Yes	FHD
FMD-190FHPL	9.2mm	2.8mm	3000mm	U180°/D130°	L100°/R100°	140°	3	Yes	FHD
FMD-1B0FHPL	11.5mm	3.2mm	3500mm	U180°/D130°	L160°/R160°	140°	3	Yes	FHD

Light Source L-10

Model	Power Consumption	Light	Air Pump Pressure	Diameter	Bright	Air Pump
FMD-10	140W	LED	45Kpa-65Kpa	∅ 8.5mm	9 Levels	3 Levels