




Veteriner Disposable Mesh Snare

Usage : It is used to grasp and remove foreign body from the alimentary tract in human body in combination with endoscope.

Special features

- Integral rotation handle
- With good flexibility and net design, it can accurately and firmly grasp the foreign body.
- Singleuse, avoiding cross contamination
- ESD and EMR's indispensable partner. The clots, polyps, tissues can be thoroughly retrieved.



The design of the handle is rotatable, the net shape grasper design avoids the loss of samples, deeper elastic nylon net, which will efficiently filter out the natural tract fluid, it is good for grasp and control.

Model	Diameter of the outer tube (mm)	Length (mm)	Actual Opening width
FMD-YWQ-TY-18x1600x20	1.8	1600	15
FMD-YWQ-TY-23x1600x20	2.3	1600	15
FMD-YWQ-TY-23x1600x30	2.3	1600	25
FMD-YWQ-TY-23x1600x40	2.3	1600	35
FMD-YWQ-TY-23x1600x20	2.3	2300	15
FMD-YWQ-TY-23x1600x30	2.3	2300	25
FMD-YWQ-TY-23x1600x40	2.3	2300	35

