



Veterinary Disposable Biopsy Forceps



Quality and ergonomics perfectly combined.

The disposable biopsy forceps from stand out due to their high quality and ergonomic design. These features enable an easy and safe removal of tissue samples or polyps. Thanks to the first class manufacture, the biopsy forceps can be inserted very smoothly into the working channel and can be controlled precisely even with sharply angled endoscopes. The extra sharp jaws can be closed accurately and thus ensure a clean cut with minimal effort.


Special features

- High-quality materials and manufacturing
- Simple and precise to use
- Sharp cutting edge for diagnostically conclusive biopsies
- Full closure of the cutting edges
- Special scissor design conserves the working channel
- Extensive product range

The extensive range offers the perfect solution for each application: whether with spike or without, with smooth or serrated jaws, with coating or without.

Ergonomic handle for optimum control.

For precise control of the jaw each biopsy forceps possess an ergonomically designed handle. It is particularly comfortable to hold, what makes the work easy and relaxed. Every the smallest rotation or movement of the handle is passed exactly to the jaw part so that this is controlled absolutely reliable.



Model	O.D. (mm)	Length (cm)	Characteristics	Endoscope Channel	Description
FMD-FB-A-S-18-2300	1.8	160-230	Oval Cup, without needle	2.0	Disposable Biopsy Forceps
FMD-FB-C-S-18-2300	1.8	160-230	Alligator Without Needle	2.0	Disposable Biopsy Forceps
FMD-FB-A-N-10-1200	1.0	1200	Oval Cup, without needle		Disposable Biopsy Forceps
FMD-FB-C-N-10-1200	1.0	1200	Alligator Without Needle		Disposable Biopsy Forceps

