

Veterinary Disposable Stone Removal Basket

ENDOcept

Special features

- Innovative handle design with push, pull and rotate functions, gallstones and foreign matter are easier to catch,
- of contrast media with the injection port on the handle. suitable for injection
- It is made of advanced alloy material and is difficult to remove. retains its shape even after stones

Purpose of usage

Gallstones from the bile duct removal and upper and lower digestion removal of foreign materials from areas



Model	O.D. (mm)	Length (mm)	Endoscope Channel (mm)		Description (mm)
FMD-RB-E4NN-1840-18120	1.0	120	1.0	Raw material: Nitinol	Disposable Basket



Model	O.D. (mm)	Length (mm)	Endoscope Channel (mm)	Characteristics	Açıklamala (mm)
FMD-RB-E4NN-1840-18120	1.8	160-230	2.0	Raw material: Nitinol	Disposable Basket



www.famed.com.tr