




FAMED SAĞLIK TIBBİ CİHAZLAR
İTH. İHRC. SAN. ve TİC. LTD. ŞTİ.



Portable USB Veterinary Endoscopy




www.famed.com.tr

info@famed.com.tr




CMOS+LED Technology

Non-fiber design makes it durable and easy for maintenance.



USB Output

Features




HD Image

720P High Definition



CMOS+LED

No Fiber



Big Channel


2.8mm




Multi-function Buttons

①White Balance ②Brightness ③Freeze

Application



Companion animals upper and lower respiratory tract.



Large animal upper and lower respiratory tract.

Specification

Model	Insertion Tube	Working Channel	Working Length	Bending	Rotation	View Angle
	4.8mm	2.0mm	1000mm	U180°, D130°	L120°, R120°	120°
	5.8mm	2.8mm	1500mm	U180°, D130°	L120°, R120°	120°