




# Disposable Video Endoskopi System



## Disposable Digital Endoscopy

Model No.	Insertion Part(mm)	Working Channel (mm)	Field of View	Depth of Field(mm)	Processor	Working Length (mm)	Deflection	Direction of View
FMD-SZQA-2812	Φ2.8	1.2	120°	5-50	SP-TXCLQ 12.1	1000	Up 180° Down 180°	0° (Forward Viewing)
FMD-SZQA-3816	Φ3.8	1.6						
FMD-SZQA-4222	Φ4.2	2.2						
FMD-SZQA-4826	Φ4.8	2.6						
FMD-SZQA-5832	Φ5.8	3.2						
FMD-SZQB-4816	Φ4.8	1.6	140°	1-80	Output (USB)			
FMD-SZQB-5826	Φ5.8	2.6						

**Deflection Angle**  
180° ± 10°



### Product Application

Diagnosis and surgery of the respiratory tract:  
Procedures such as lesion observation and biopsy



# Disposable Video Endoskopi System



- Photo and video recording with an external display
- Electronic imaging with a clearer and more realistic image
- Avoid cross-infection and improves surgical safety
- Avoid poor observation due to product aging and performance degradation
- Reducing the cost of use, no need for disinfection and maintenance for the insertion tube
- Significantly shorten the procedure and articulation time to improve efficiency

## Processor

Processor Model	SP-TXCLQ12.1
Power Supply Input	DC 12V
Screen Size	12.1"
Resolution	1280*800
Storage	≥16G
External Interface	Type-C, HDMI Output, three-phase power interface
Screen Type	Capacitive touch screen



Processor Model	SP-TXCLQ-02-3.5
Design	Portable
Screen Type	3.5" IPS Screen
Screen Rotation Angle	Up160° , down 30°
Resolution	640×480
Battery	2000mah
External Interface	Type-C, MINI HDMI



**FAMED SAĞLIK TIBBİ CİHAZLAR İTH. İHRC. SAN. ve TİC. LTD. ŞTİ.**

**Address::** Cevizli Mh. Hacılar Cd. No:44/A Kartal - İstanbul / TURKEY

**Phone :** +90 216 597 0 567 / +90 216 441 29 99 **Fax:** +90 216 397 15 15

**W:** www.famed.com.tr **E:** info@famed.com.tr

| [www.famed.com.tr](http://www.famed.com.tr)