



BIORAD®

ENZYM PERFECT

Pre-cleaning And
Disinfection



FORMS OF PACKAGE

Description	Package
5 L PE drum	4 x 5 L
10 L PE drum	1 pc.

CHEMICAL COMPOSITION:

5-enzyme ((Amylase, Lipase, Protease, Cellulase, Mannanase), Nonionic tensides, Benzalkonium chloride, Corrosion inhibitors and excipients.

PRODUCT SPECIFICATIONS:

ENZYM PERFECT is a ready-to-use instrument cleaner and disinfectant which contains enzyme. It is used in the enzymatic cleaning and disinfection of medical and surgical instruments and heatresistant and heat-sensitive materials. It cleans the organic dirt on the instruments and disinfects at the same time. It prevents corrosion. It does not contain aldehyde, phenol and phosphate.

AREAS OF USE:

It is suitable for use for enzymatic cleaning and disinfection of endoscopes, anesthetic apparatus, dental instruments, heat-resistant and heat-sensitive materials, instrument disinfection tubs, manual, ultrasonic washing machines and automatic washing machines including treated medical and surgical instruments and MIC surgical instruments as well as all the tools made of ceramic, porcelain and glass. It is recommended to use as excipient in the cleaning and disinfection of Wolf, Storz, Fujinon, Pentax, Olympus branded endoscopes.

USE:

Sufficient amount of solution is used without dilution. The product is ready-to-use

ANTI-MICROBIAL CHARACTERISTICS:

Bactericide, Fungicide, Virucide, Mycobactericidal