



BIORAD[®]

STRONG A

Washing And
Disinfection Machines



FORMS OF PACKAGE

Description	Package
5 L PE drum	4x5 L

CHEMICAL COMPOSITION:

Phosphoric acid, Benzalkonium chloride, Anti-corrosive agents

PRODUCT SPECIFICATIONS:

It is a neutralizing solution which contains liquid and organic acid and used for all washerdisinfectors. It is an effective product in the neutralization of remaining alkaline residues and removal of the potential lime scales inside the washing machines in main washing procedures. pH=1-3.

AREAS OF USE:

Strong A is used in the neutralization of washing machines in thermal and chemothermal process after the alkaline cleaning. It contains organic acid. It is used for the neutralization (40-60°C) of the remaining alkaline residues and removal of lime scales after the alkaline cleaning of the surgical instruments, rigid endoscopes, MIC instruments, anesthetic instruments in the hospitals and medical institutions.

USE:

BIORAD STRONG A is suitable for use in all washer-disinfectors. The recommended usage concentration is 0.1-0.3% depending on the rate of contamination and water quality. Follow the recommendations of the machine manufacturer for dosing information