



BIORAD®

CLEAN PERFECT

Surface and Floor
SPRAY
Disinfectants



FORMS OF PACKAGE

Açıklama	Ambalaj
1 L PE bottle	20×1 L
5 L PE drum	4×5 L

CHEMICAL COMPOSITION:

Didecyldimethylammoniumchloride , Ethanol, Corrosion inhibitor and excipients

PRODUCT SPECIFICATIONS:

It is an alcohol based, anti-corrosive product which has a pleasant odor and used for the quick disinfection of medical device surfaces. It is dermatologically tested. It is used as an aid in the quick disinfection of all medical devices, manual instruments, equipment and their surfaces which are alcohol-resistant and suitable for moistening.

AREAS OF USE:

It is used for the cleaning and disinfection of all horizontal and vertical surfaces and especially PVC epoxy areas, floors, walls, ultrasonography probes, incubators, aerators and micromotor tips used in dentistry, stethoscopes, cabinets, doors, examination, disinfection requiring medical devices and equipment as well as operating rooms, intense care rooms, sickbed edges, iv holders, areas contacting with hands and skin on microscopes and tools contacting with patients' skin in the hospitals that require quick disinfection.

USE:

Clean Perfect Spray, is a ready-to-use spray disinfectant suitable for use on all horizontal and vertical surfaces. It is sprayed from a distance of 30cm-40cm . It takes 30 seconds for the effect to appear. It is wiped with a clean cloth and allowed to dry.

ANTI-MICROBIAL CHARACTERISTICS:

Bactericide, Fungicide, Virucide, Antimycobacterial, Tuberculocide

www.famed.com.tr