



# BIORAD®

## READY SPRAY WIPE

Surface And Floor  
Disinfectants



### FORMS OF PACKAGE

Description	Package
25 Piece	

#### CHEMICAL COMPOSITION:

Benzalkonium chloride, Didecylmethylammoniumchloride and excipients

#### PRODUCT SPECIFICATIONS:

Ready Spray Wipe is ready-to-use and has an extended antimicrobial action spectrum with the quaternary ammonium compounds in its composition. It does not contain aldehyde, phosphate and alcohol. It is recommended to use on surfaces with metal parts since it contains corrosion inhibitor. Ready Spray Wipe is compatible with glass, ceramic, silicone, plastic (including Plexiglas), wood and stainless steel materials.

#### AREAS OF USE:

It is used for the cleaning and disinfection of all horizontal and vertical surfaces and especially incubators, PVC epoxy areas, floors, walls, ultrasonography probes, incubators, aerators and micromotor tips used in dentistry, disinfection requiring medical devices and equipment in ENT units as well as operating rooms, intense care rooms, in the hospitals that require quick disinfection

#### ANTI-MICROBIAL CHARACTERISTICS:

Bactericide, Fungicide, Virucide, Antimycobacterial, Tuberculocide