



BIORAD® AQUA OXYGEN

Endoscope And Medical
Instrument Disinfection



FORMS OF PACKAGE

Description	Package
1 litre PE bottle	20 x 1 litre
5 litre PE drum	4 x 5 litre
10 litre PE drum	1 piece

CHEMICAL COMPOSITION:

Hydrogen peroxide deionised water

PRODUCT SPECIFICATIONS:

Aqua Oxygen a combined, hydrogen-peroxide based disinfectant for disinfection of the medial equipment and it is aldehyde-free and phenol-free. It contains corrosion inhibitors.

INTENDED USE:

Aqua Oxygen is used for cleaning and disinfection of surgical equipment and septic scots and root canals in the dentistry. It is used for surgical purpose and cleaning of scots and root canals in the dentistry as antiseptic-disinfectant and for removal of ear wax as ear drop. It is compliant with glass, ceramic, silicone, plastic, acrylic, polymer, polycarbonate, polystyrene, HDPE, porcelain and stainless steel

ANTI-MICROBIAL CHARACTERISTICS:

Bactericidal, Fungicidal, Virucidal, Tuberculocidal, Micobactericidal