




BIORAD®

HIPO

Endoscope And Medical
Instrument Disinfection



FORMS OF PACKAGE

Description	Package
1 litre PE bottle	20 x 1 litre
5 litre PE drum	4 x 5 litre
10 litre PE drum	1 piece

CHEMICAL COMPOSITION:

Sodium hypochloride

PRODUCT SPECIFICATIONS:

It acts and kills enzymes which play key role in the chlorine cell

INTENDED USE:

Hypo is used for cleaning and disinfection of medical and surgical equipment. It is compliant with metals, plastics, rubber, glass, wood, porcelain, glass, ceramic, silicone, plastic, acrylic, polymer, polycarbonate, polystyrene, HDPE porcelain and stainless steel.

USAGE CONCENTRATION

It is diluted by 1/10 or 1:100 depending on amount of organic substance on the surface to be disinfected. Dilution of 1/10 is recommended for disinfection of blood that fell on the floor or splashed around. 5ml is added per 1L for washing machine or manual washing. Its short time action is 5 minutes.

ANTI-MICROBIAL CHARACTERISTICS:

Bactericidal, Fungicidal, Virucidal, Tuberculocidal, Micobactericidal