

# HEMATOLOGY REAGENTS



## FAMED CLEANER M-30P REACTIVE CLEANING CHEMICALS

FAMED CLEANER M-30P Reactive Cleaning Chemicals (Rinse, Cleaner) It is a type of detergent used for automatic cleaning of internal hoses (tubing), baths and other necessary parts of the device. It also functions as a means of keeping the hydraulic system of the device wet. The device automatically cleans itself after each sample count and becomes ready for the next count. 8 Enzymatic cleaners are enzyme-based, isotonic cleaning solutions. It is used for cleaning and wetting tubing and mixing baths. It is prepared in a way that will not damage the properties of the plastic parts of the device it comes into contact with.



### BRAND HEMATOLOGY REAGENTS

Mindray Serisi- Medica Serisi- Konelab Serisi - Erba Serisi – Beckman Serisi – Dirui Serisi – Olympus Serisi – Abbott Serisi – Roche Serisi – Randox Serisi – Human Serisi – Siemens Serisi – Medonic Serisi – Swelab Serisi – Hemacell Serisi – Abx Micros Serisi – Diatron serisi -Sysmex -Pentra serisiNihon Kohden Celltac Serisi - Sysmex 3-diff Serisi- Horiba ABX Micros 60/45 SerisiSysmex XS/ XT/ XE Serisi-

**FAMED SAĞLIK TIBBİ CİHAZLAR İTH. İHRC. SAN. ve TİC. LTD. ŞTİ.**

**Address::** Cevizli Mh. Hacılar Cd. No:44/A Kartal - İstanbul / TURKEY

**Phone :** +90 216 597 0 567 / +90 216 441 29 99 **Fax:** +90 216 397 15 15

**W:** www.famed.com.tr **E:** info@famed.com.tr

✉ | [info@famed.com.tr](mailto:info@famed.com.tr)

🔍 | [www.famed.com.tr](http://www.famed.com.tr)