



FAMED SAĞLIK TIBBİ CİHAZLAR  
İTH. İHRC. SAN. ve TİC. LTD. ŞTİ.


## Özellikleri

- 10 x Büyütme
- 32 mm geniş görüş alanı
- Polarize, polarize olmayan ve amber ışık aydınlatma
- Çıkarılabilir koruyucu cam
- Otomatik kapanma
- Birlikte verilen adaptör tüm telefonlara uyar
- Tamamı metal gövde

## DERMATOSKOP DE-3100



## DERMATOSKOP DE-3100



### Teknik Bilgiler

<b>Malzeme</b>	Alüminyum ve Cam
<b>Optik Tasarım</b>	Tamamı cam, 4 elemanlı 3 grup, yansıma önleyici kaplama
<b>Mercek Çapı</b>	32 mm (ön); 25mm(arka)
<b>Büyütme</b>	10X
<b>Çarpıtma</b>	8%
<b>Polarizasyon</b>	4 LED beyaz polarizesiz, 4 LED sarı polarize, 8 LED beyaz polarize
<b>Çözünürlük</b>	300 LP/MM (Eksen) 250 LP/MM (Kenar)
<b>Led Tipi</b>	SMD LED boncuklar
<b>Pil Kapasitesi</b>	1000mAh Lityum iyon
<b>Dolum</b>	USB-C
<b>Çalışma Zamanı</b>	3-6 Saat



[www.famed.com.tr](http://www.famed.com.tr)



[info@famed.com.tr](mailto:info@famed.com.tr)



## Accessories

# DE-3100



### UNIVERSAL PHONE ADAPTER

Save images of lesions on your smartphone



### IN THE BOX

Relative accessories included in the package



### SMALL CONTACT PLATE

To detect difficulty access pigmented lesions



### CAMERA ADAPTER

Save images of lesions on your camera



### PHONE COVER ADAPTER

iPhone & Samsung phone cover for dermascope



### BACK UP BATTERY

1000mAh Lithium ion battery for replacement

