



FAMED SAĞLIK TIBBİ CİHAZLAR
İTH. İHRC. SAN. ve TİC. LTD. ŞTİ.


Features

- 10 x Magnification
- 32mm wide field of view
- Polarized, non-polarized, and amber light illumination
- Detachable protective glass
- Automatic shutdown
- Included adapter fits all phone
- All-metal housing

DERMATOSCOPE DE-3100



DERMATOSCOPE DE-3100



Technical Information

Material	Aluminum & Glass
Optical Design	All glass, 4 elements 3 groups, anti-reflection coating
Lens Diameter	32mm(front); 25mm(rear)
Magnification	10X
Distortion	8%
Polarization	4 LEDs white non polarized, 4 LEDs amber polarized, 8 LEDs white polarized
Resolution	300 LP/MM (Axis) 250 LP/MM (Edge)
LED Type	SMD LED beads
Battery Capacity	1000mAh Lithium ion
Charging	USB-C
Working Time	3-6 hours



www.famed.com.tr



info@famed.com.tr



Accessories

DE-3100



UNIVERSAL PHONE ADAPTER

Save images of lesions on your smartphone



IN THE BOX

Relative accessories included in the package



SMALL CONTACT PLATE

To detect difficulty access pigmented lesions



CAMERA ADDAPTER

Save images of lesions on your camera



PHONE COVER ADAPTER

iPhone & Samsung phone cover for dermascope



BACK UP BATTERY

1000mAh Lithium ion battery for replacement



Address:: Cevizli Mh. Hacilar Cd. No:44/A Kartal - İstanbul / TURKEY
Phone : +90 216 597 0 567 / +90 216 441 29 99 **Fax:** +90 216 397 15 15
W: www.famed.com.tr **E:** info@famed.com.tr