



# DERMATOSCOPE DE-4100



FAMED SAĞLIK TIBBİ CİHAZLAR  
İTH. İHRC. SAN. ve TİC. LTD. ŞTİ.

## Technical Information

Material	Aluminum & Glass
Optical Design	All glass, 4 elements 3 groups, anti-reflection coating
Lens Diameter	33mm(front); 33mm(rear)
Magnification	10X
Distortion	8%
Polarization	6 LEDs white non polarized, 4 LEDs amber polarized, 12 LEDs white polarized
Resolution	300 LP/MM (Axis) 250 LP/MM (Edge)
LED Type	SMD LED beads
Battery Capacity	1000mAh Lithium ion
Charging	USB-C
Working Time	3-6 hours



## Features

- 10 x Magnification
- 32mm wide field of view
- Polarized, non-polarized, and amber light illumination
- Detachable protective glass
- Automatic shutdown
- Included adapter fits all phone
- Naked eye detection
- Bright LED illumination lights with 3 levels of brightness control
- All-metal housing

[www.famed.com.tr](http://www.famed.com.tr)

[info@famed.com.tr](mailto:info@famed.com.tr)



# Accessories DE-4100



**UNIVERSAL PHONE ADAPTER**  
Save images of lesions on your smartphone



**IN THE BOX**  
Relative accessories included in the package




**SMALL CONTACT PLATE**  
To detect difficulty access pigmented lesions



**CAMERA ADAPTER**  
Save images of lesions on your camera



**PHONE COVER ADAPTER**  
iPhone & Samsung phone cover for dermascope



**BACK UP BATTERY**  
1000mAh Lithium ion battery for replacement



**Address::** Cevizli Mh. Hacilar Cd. No:44/A Kartal - İstanbul / TURKEY  
**Phone :** +90 216 597 0 567 / +90 216 441 29 99 **Fax:** +90 216 397 15 15  
**W:** www.famed.com.tr **E:** info@famed.com.tr