

# DERMATOSKOP PRO DE-4100




FAMED SAĞLIK TIBBİ CİHAZLAR  
İTH. İHRC. SAN. ve TİC. LTD. ŞTİ.


## DERMATOSKOP PRO

### Technical Information

|                  |   |
|------------------|---|
| Material         | Aluminum & Glass  |
| Optical Design   | All glass, 4 elements 3 groups, anti-reflection coating         |
| Lens Diameter    | 33mm(front); 33mm(rear)   |
| Magnification    | 10X   |
| Distortion       | 8%  |
| Polarization     | Polarization & non-polarization & UV<br>12 LEDs white polarized |
| Resolution       | 300 LP/MM (Axis) 250 LP/MM (Edge)                               |
| LED Type         | SMD LED beads   |
| Battery Capacity | 1000mAh Lithium ion   |
| Charging         | USB-C   |
| Working Time     | 3-6 hours   |



Polarizasyonsuz



UV

### Features

- 10 x Magnification
- 32mm wide field of view
- Polarized, non-polarized, and amber light illumination
- Detachable protective glass
- Automatic shutdown
- Included adapter fits all phone
- Naked eye detection
- Bright LED illumination lights with 3 levels of brightness control
- All-metal housing



[www.famed.com.tr](http://www.famed.com.tr)



[info@famed.com.tr](mailto:info@famed.com.tr)



## DERMATOSKOP PRO Accessories DE-4100



**UNIVERSAL PHONE ADAPTER**  
Save images of lesions on your smartphone



**IN THE BOX**  
Relative accessories included in the package




**SMALL CONTACT PLATE**  
To detect difficulty access pigmented lesions



**CAMERA ADDAPTER**  
Save images of lesions on your camera



**PHONE COVER ADAPTER**  
iPhone & Samsung phone cover for dermascope



**BACK UP BATTERY**  
1000mAh Lithium ion battery for replacement



**Address::** Cevizli Mh. Hacilar Cd. No:44/A Kartal - İstanbul / TURKEY  
**Phone :** +90 216 597 0 567 / +90 216 441 29 99 **Fax:** +90 216 397 15 15  
**W:** [www.famed.com.tr](http://www.famed.com.tr) **E:** [info@famed.com.tr](mailto:info@famed.com.tr)