




FAMED SAĞLIK TIBBİ CİHAZLAR
İTH. İHRC. SAN. ve TİC. LTD. ŞTİ.

WOODS LAMP DE-215



UV Light



White Light



What It Has

- **365nm UV Light:** A Wood's lamp has broadband light sources that emit light at wavelengths between 320nm and 400nm, with a peak at 365nm.
- **4.5X Magnification:** Low distortion magnify
- **2000mAh Battery:** Offering up to 6 hours long time and stable diagnosis.

Specs

Lens Diameter	50mm
Magnification	4.5X
Ultraviolet Wavelength	365nm
Radiation Intensity	3.5 mW/cm ²
Battery Capacity	2000mAh
Charging	USB Type-C
Working Time	2-6 hours
Dimension	240mm*100mm*34mm(L*W*H)
Weight	166.5g