




# Disposable Stone Extraction Basket

Everything in a safe grip




It is used to grasp and remove foreign body from alimentary tract in human body in the combined use of endoscope.

## Special features

- Integral rotation handle
  - With good flexibility and net design, it can accurately and firmly grasp the foreign body.
  - Disposable use, avoiding cross contamination
- ESD and EMR's indispensable partner. The clots, polyps, tissues can be thoroughly retrieved.

## Efficient control by 3-ring handle.

The 3-ring handle combines several advantages. Thanks to its ergonomic design the handle is extremely user-friendly and easy to operate with one hand. At the same time it relays every movement and rotation exactly to the forceps, grasper or basket, what makes the intervention even safer.



Model	Diameter of the outer tube (mm)	Length (mm)	Opening width
FMD-YWQ-ZJ-18x1600x15	1.8	1600	15
FMD-YWQ-ZJ-23x1600x15	2.3	1600	15
FMD-YWQ-ZJ-23x1600x25	2.3	1600	25
FMD-YWQ-ZJ-23x1600x35	2.3	1600	35
FMD-YWQ-ZJ-23x2300x15	2.3	1600	15
FMD-YWQ-ZJ-23x2300x25	2.3	1600	25
FMD-YWQ-ZJ-23x2300x35	2.3	1600	35

