





EndoPik

Laparoscopic Video Laparoscopy System

VES-100



Characteristics

- All in one system (CCU + Light source + Recorder)
- Saving cost and efficient use based on " All in one system "
- Power LED light source brighter than Xenon 300W light source
- Saving and recording images into SD card or HDD
- Touch LCD and Intuitive GUI for users



"ALL IN ONE SYSTEM!"



SOMETECH
SINCE 1989




SAĞLIK TIBBİ CİHAZLAR ve SAN. TİC. LTD. ŞTİ.



Odak ▲/▼ Zoom IN/OUT.
Capture button



One touch
coupler




Camera

System Specifications

Category	Specification
Size (WxDxV)	310mm x 300mm x 150mm
Voltage / Frequency	AC 100 - 240V, 50 / 60 Hz
Video output	HDMI
Video source	Color, Brightness, Contrast, Sharpness, Exposure area, Exposure level.
White balance	Push set
Zoom	Digital Zoom (1.0-2.0 X)
Image rotation	180°
GUI interface	TFT Color LCD
Light cable length	3.5M
Light Source	White LED (Cool light)
Light intensity	115 mm'de 80,000 Lux
Light adjustment	0(OFF) 1 - 10 (Level Adjustment, 10= Max)
Life Time	60.000 saat
Color Temperature	6800 K
Still image capture	Full HD (1920X1080) .jpg
Video recording	MP4(H.264), .avi
Storage media	SDHC/USB2.0

- Crystal clear FULL HD(1080P) resolution
- Durable and convenient button on camera head
- Easy zooming and fine focusing by button on camera head
- Compact, convenient and light-weight camera head
- Water proof camera head



Touch LCD and user - Friendly GUI

Still and video images in SD card



Image saving to HDD touch USB part

Camera

Category	Specification
Resolution	FULL HD 1920 X1080P
Sensor	1/3" CMOS sensörü
Camera & Coupler	Waterproof camera head and C-Mound coupler
Sterilization	Low temperature sterilization (E.O Gas, Plasma, Chemicals)
Cable length	3.5M

www.famed.com.tr



Cevizli Mah. Hacılar Cad.
No:44/A Kartal / İstanbul



+90 216 597 0 567
+90 216 441 29 99



+90 216 397 15 15